


CE 0483

Triceram®

on zirconium oxide

The ceramic
for CAD/CAM technology

Product information and instructions for use



Dear customer,

your choice of Dentaaurum products welcomes you to a world of well harmonized and meticulously coordinated quality products for the production of aesthetic ceramic work.

The most important factor when producing ceramic restorations is to maintain precision whilst following the Triceram® instructions for use for zirconium oxide.

In this brochure you will find many practical tips for using our products. Further information for processing Triceram® can be found in our brochure "Working tips for Triceram®" (REF 989-676-20).


Additionally, we are here to support you should you have any more questions on working with our products.

Our dental technical advisors are available for any questions regarding working with our products.

Hotline Tel. No. +497231/803-440

Information and instructions for use on the Dentaaurum ceramic system can be found in the internet under

www.dentaaurum.de



Contents

1. Triceram®, the ceramic for CAD/CAM technology

1.1. Zirconium oxide (ZrO ₂)	4
1.2. Bonding strength, flexural strength and chemical resistance	4
1.3. Tips for the ceramic bond	5

2. Processing techniques

2.1. Framework design and preparation	6
2.2. Advice on processing the framework	7
2.3. Applying the Fluo Liners	7
2.4. Using the shoulder material	8
2.5. Applying the shoulder material	8
2.6. Layering example	9
2.7. Layering step-by-step	10
2.8. Using the Stains Universal	12

3. Information about the ceramic system

3.1. Products and their components	13
3.2. Shade classification table	14
3.3. Firing test	14
3.4. Cleaning the furnace	15
3.5. Firing table standard program	16
3.6. Firing tables for selected furnaces	17

4. Individual products 20

5. Symbol description 23

1. Triceram®, the ceramic for CAD/CAM technology

1.1. Zirconium oxide (ZrO₂)

Most of the zirconium oxide (ZrO₂) used in restorative dentistry is so-called yttrium stabilized TZP ceramic (tetragonal zirconia polycrystal). This is a polycrystalline ceramic consisting of yttrium stabilized tetragonal ZrO₂ particles. This material decreases outstandingly high mechanical strength and fracture toughness, it is white or off-white in colour with a certain translucency resembling natural dentition and opens up entirely new possibilities for large spanned, metal-free, aesthetical and challenging restorations. The high strength and fracture toughness of zirconium oxide ceramic (Y-TZP) is estimated to withstand extreme, long-term stress in the molar region with excellent clinical prognosis.

The outstanding strength of TZP ceramics is due to its fine crystalline structure and ability to resist cracking through transformation. Y-TZP ceramics consist of tetragonal zirconium oxide particles, which due to their size, the pressure of the surrounding matrix and the addition of Yttrium oxide, remain stable at room temperature. When these tetragonal stabilized particles are subjected to energy in the form of tension, such as from a progressing crack, they transform into the stable monocline modification. The material is able to absorb the tension and due to a growth in volume, which goes in hand with the tetragonal monocline transformation, compressive strength is increased in that particular area. Not only is the tension removed from the crack tip but the build up in compressive strength increases the resistance.

1.2. Bonding strength, flexural strength and chemical resistance

Triceram® is the ideal partner for the materials titanium and zirconium oxide in the dental CAD/CAM technology. Both framework materials produce an exceedingly strong bond by means of special procedures. Zirconium oxide achieves a bonding strength of 40 MPa (according to Schmitz-Schulmeyer).




Fig. 1



Fig. 2



Fig. 3



Fig. 4

The fluo liner which has been specially developed for zirconium oxide guarantees on white, unstained framework the accurate reproduction of vita classic tooth shades and natural fluorescence without the hindrance of light transmission from the framework. The special physical properties of Triceram® on zirconium oxide framework and its favourable aesthetics have been confirmed in various independent studies. In addition Triceram® clearly beats the required minimal standard (EN ISO 6872) for flexural strength (Fig. 5) and considerably exceeds the required level of chemical resistance (Fig. 6).




Fig. 5: Flexural strength

Solubility

The permissible level of solubility in the hydrolysis test is max. 100 µg/cm² according to ISO 9693.

Measured values Triceram®:

Opaque: 55 µg/cm²

Dentine: 31 µg/cm²

Incisal: 31 µg/cm²

This shows the solubility is distinctly lower than the permissible level.

Fig. 6: Solubility of individual materials

CTE (25-500 °C / 77-932 °F):

Triceram® Liner: → 9.4 x 10⁻⁶ K⁻¹

Triceram® Dentine: → 8.6 x 10⁻⁶ K⁻¹

Triceram® Incisal: → 8.6 x 10⁻⁶ K⁻¹

Fig. 7: CTE values of individual materials

1.3. Tips for the ceramic bond

The physical and thermal properties of Y-TZP and its unique appearance necessitates a specially developed veneering ceramic and processing procedure. As with other familiar ceramic materials, the zirconium oxide framework and its relevant ceramic forms a bonds primarily using compressive strength. For this reason the CTE of the ceramic Triceram® (Fig. 7), must be smaller than the CTE of the zirconium oxide framework (CTE approx. 10.5 x 10⁻⁶ K⁻¹). In order for the ceramic to shrink evenly and develop an equal amount of compressive strength across the substructure, the framework must be designed as a smaller anatomical version of the end product, without any sharp sides or edges. The veneering ceramic requires sufficient space in order to ensure it is applied in an even thickness. Thicknesses of over 2 mm must be avoided. Please observe the information in section "processing techniques" in these instructions for use.

- ¹
 - University of Freiburg, Prof. Dr. Kappert H.F.: Bruchfestigkeit und Frakturverhalten von DCS-gefrästen Frontzahnkronen nach Verblendung mit Vitadur D und Triceram®, ZWR, 2001: 134-139, (DC-Zirkon®)
 - University of Nantes, Frankreich, Praud C.: Dissertation: Rapport de la Céramique Triceram® à la réalisation de prothèses céramo-métalliques sur Titane: 1999
 - University of München, Bader R.-D., Aschl I., Rammelsberg P. Pospiech P.: IADR Abstracts 1021: Improved bond strength of low fusing ceramics to titanium: 2000
 - University of Zürich, Dr. Irena Sailer, Dr. Claudia Holderegger, Dr. Roland Jung, ZTM Ana Suter: Zirkonoxid-Verblendkeramiken: Farbstabilität und technische Verarbeitung, Quintessenz Zahntechnik 5/2005, S. 498-512
 - University of Modena, ZTM Bruno Marziali, Dr. Giuliano Malaguti und Prof. Sergio Bortolini: Titan und Zirkoniumdioxid, Quintessenz Zahntechnik 6/2006, S. 632-647
 - University of Zürich, PD Dr. med. dent. Dr. rer. nat. Jens Fischer, Dipl.-Ing. (FH) Bogna Stawarczyk: Zirkoniumdioxid und Titan: Zwei Gerüstwerkstoffe – eine Verblendkeramik, DZW-Zahntechnik 12/06, S. 20-22

2. Processing techniques

2.1. Framework design and preparation

The dental ceramic and framework materials are usually designed to harmonize thermally. This ensures that the state of stress, which occurs during the cooling phase after having been fired, is ideal for both materials. Amongst other factors, the state of stress is influenced by the geometrical design of the substructure.

For this reason, please observe the following basic principles:

- The framework must represent the final tooth shape in reduced anatomical form. Missing material must be replaced by building up the zirconium oxide framework and not the veneering ceramic.
- Avoid sharp sides and edges. Always integrate anatomically rounded forms and smooth transitions into the design. Avoid undercuts.
- The framework should be designed so that the ceramic may be applied in a uniform thickness.
- The maximum ceramic thickness is 2 mm.

Frameworks with holes or cracks must not be used!

2.2. Advice on processing the framework

Finishing and sandblasting milled or ground framework carries the risk of superficial alteration in the structure of the Y-TZP and should be reduced to a minimum. Please only use the tools recommended (eg. the finishing set, REF 137-840-00) at the relevant revolutions per minute, under moderate pressure. Take exceptional care when grinding the connectors. Please bear in mind the minimum cross section required, stipulated by the material manufacturers. When separating the fired ceramic interdentally with a separating disc, never allow the disc to cut into the underlying framework. Please do not overheat the framework when grinding, use water to cool the object down when necessary. When sandblasting the framework, please follow the manufacturers instructions. The framework should be cleaned before the ceramic build-up using a steam cleaner.




Fig. 1




Fig. 2



Fig. 3: Finishing set REF 137-840-00




Fig. 1



Fig. 2



Fig. 3



Fig. 4

2.3. Applying the Fluo Liners

Fig. 1:

Prepare the framework surface according to the manufacturers instructions.

Fig. 2:

Mix the liner to the relative tooth shade (see the mixing recommendations) using Universal B.O.L liquid REF 299-180-40, until it has a creamy consistency. Apply the Fluo Liner evenly onto the framework enough to cover it thoroughly and then fire. The Fluo Liner gives the framework the ideal base shade and fluorescence. A second coating may be given.

Mixing recommendations:

Tooth shade	FL 1	FL 2	FL 3	FL 4	FL 5	FL 6
A1	1/3				2/3	
A2	2/3				1/3	
A3	1					
A3.5	3/4					1/4
A4	1/2					1/2
B1		1/3			2/3	
B2		2/3			1/3	
B3		3/4			1/4	
B4		1				
C1			1/3		2/3	
C2			2/3		1/3	
C3			3/4		1/4	
C4			1			
D2				2/3	1/3	
D3				3/4	1/4	
D4		1/3		2/3		

Fig. 3:

After firing, the liner has a silky, shiny surface.

Fig. 4:

The liner can also be used on stained zirconium oxide framework to achieve an optimal shade and natural fluorescence.

2.4. Using the shoulder material

Triceram® offers four shoulder materials subdivided into the shade groups A, B, C and D. Using the shoulder material „white“, all shade nuances from A1 to D4 can be mixed according to the mixing table. The shoulder material „transparent“ increases the transparency within the shoulder. Further individual modifications can be achieved using the shoulder modifications “yellow” and “orange”. Use the universal shoulder material liquid for mixing.

Mixing recommendations:

Tooth shade	A	B	C	D	white
A1	1/2	-	-	-	1/2
A2	2/3	-	-	-	1/3
A3	3/4	-	-	-	1/4
A3.5	1/1	-	-	-	-
A4	1/1	-	-	-	-
B1	-	1/2	-	-	1/2
B2	-	2/3	-	-	1/3
B3	-	1/1	-	-	-
B4	-	1/1	-	-	-
C1	-	-	2/3	-	1/3
C2	-	-	3/4	-	1/4
C3	-	-	3/4	-	1/4
C4	-	-	1/1	-	-
D2	1/3	-	-	1/3	1/3
D3	1/2	-	-	1/4	1/4
D4	-	-	-	1/1	-

2.5. Applying the shoulder material

Isolate the area to be built-up with the separating agent SM-Isokit (REF 260-324-01). Apply shoulder material from the framework to the preparation line. Allow the material to dry (using dryer or tissue), remove the coping and fire using the appropriate firing program.

The shoulder material shrinks during the sinter process in the first bake. A second and third firing will correct the material reduction. Apply SM-isokit to the areas to be built-up and then re-apply the missing material. Fire exactly the same as the first bake. After the second bake has been completed fit the coping back onto the die using the recommended grinding burs (watch the rpm! Max. 15 000 U/min). Carry out the ceramic build-up in the conventional manner.



Fig. 1: Shortened coping, after the opaque bake




Fig. 2: 1. Applying the shoulder material



Fig. 3: Coping before the 2. shoulder baks, shoulder material has been built-up



Fig. 4: Coping after fitting




Fig. 1: Standard build-up



Fig. 2: Individual build-up



Fig. 3: Individual build-up



Fig. 4: Standard build-up

- FL Fluo Liner
- OD Opaque Dentine
- D Dentine
- IT Incisal transparent
- IO Incisal opal
- NT Neutral transparent
- Modifier (Fluo Dentine, Chroma Dentine)



Fig. 5: Alternative layers in the incisal area



Fig. 6: Detailed view

2.6. Layering example

Select the ceramic material according to the required tooth shade. Follow the layering suggestions in the layering diagrams.

2.7. Layering step-by-step

Fig. 1:

Build-up the reduced anatomical tooth form with opaque dentine.

Apply the opaque dentine OD over the entire labial surface. Ensure that the OD overlaps the incisal edge therefore lengthening the ZrO_2 substructure. (increases the flow of light)

Fig. 2:

Layer the dentine over the opaque dentine.

Fig. 3:

Build-up the entire anatomical shape with the dentine.

Fig. 4:

Cut back the dentine systematically in order to achieve enough space for an optimal incisal layer thickness.

Fig. 5/9:

Complete the incisal area by building up with transparency material (IT 57-IT 60).

The finished build-up before firing. With bridgework always be sure to separate the ceramic between the teeth interdentally, down to the ZrO_2 framework, before the first firing to allow for controlled shrinkage.

Fig. 6:

The ceramic crown after the first dentine bake. The fired results show a silky matt surface. The incisal edge is sharp and defined.

Fig. 7:

In the correction bake after the first dentine bake the dentine is applied so that it tapers off towards the incisal edge. This is then built-up in incisal material.

Fig. 8:

Grinding the ceramic: Grind the shape and correct the function using diamond sintered burs or stones. The surface conditioning and surface design can also be undertaken using diamond sintered burs, stones or rubber burs. After grinding, carefully steam and clean.

Fig. 10/11:

The object is glazed without the use of glaze material. Paint specific colour modifications directly onto the surface with stains. Mix the stains using stains universal liquid. Conduct the glaze bake as required. Rubbered surfaces require lower final temperature. A natural looking gloss can also be achieved by polishing the surface.




Fig. 1



Fig. 2



Fig. 3



Fig. 4




Fig. 5



Fig. 6




Fig. 7




Fig. 8




Fig. 9




Fig. 10




Fig. 11

2.8. Using the Stains Universal

The stains universal allow shade modifications to be made. They can either be embedded within or mixed together with the ceramic. It is possible to mix up to 10 % stain into various ceramic materials. In addition, these materials can be used together with the stains universal liquid to stain the surface of the ceramic. These stains universal materials are distinguished by the individual characteristics possible and depth of colour.

Body ST A	ST 7 dark pink
Body ST B	ST 8 purple
Body ST C	ST 9 blue
ST 0 neutral	ST 10 grey
ST 1 white	ST 11 olive green
ST 2 vanilla	ST 12 olive yellow
ST 3 yellow	ST 13 middle brown
ST 45 orange plus	ST 14 reddish brown
ST 6 pink	ST 15 black

Photo credits:

Bruno Martin, Le Plessis Grammoire, France

Photo sequence: The tooth and its journey through time




Fig. 1: The young tooth



Fig. 2: The middle aged tooth



Fig. 3: The mature tooth



Fig. 4: The aged tooth



Fig. 1: Master set ZR




Fig. 2: Shade wheel Basic materials I



Fig. 3: Shade wheel Effect materials II




Fig. 4: Shade wheel Effect materials III

3. Information about the ceramic system

3.1. Products and their components

FL Fluo Liner

1, 2, 3, 4, 5, 6

SM Shoulder

A, B, C, D

SM Modifier

white, transparent

OD Opaque Dentine

A1, A2, A3, A3.5, A4, B1, B2, B3, B4,
C1, C2, C3, C4, D2, D3, D4

D Dentine

A1, A2, A3, A3.5, A4, B1, B2, B3, B4,
C1, C2, C3, C4, D2, D3, D4

Bleach

Bleach CO 1, Bleach CD 1, CD 2, Bleach CI 1

CD Chroma Dentine

A, B, C

IFD Intensive Fluo Dentine

cream, yellow, orange, orange plus

IT Incisal transparent

57, 58, 59, 60

IO Incisal opalescence

57, 58, 59, 60, ice white, polar white,
snow white, sky blue, amber

NT Neutral transparent

NT

CM Correction

CM

TE Transpa Effect

white, yellow, pink, honey, blue, brown

G, GD Gingival

light, dark

ST Stains universal

A, B, C, 0, 1, 2, 3, 45, 6, 7, 8, 9,
10, 11, 12, 13, 14, 15

Universal liquid

Universal pastes, Shoulder materials,
Stains universal, Modelling liquid MV,
Modelling liquid LV standard,
Modelling liquid LV+, B.O.L.-, SM-Isosfit

3.2. Shade classification table

	A1	A2	A3	A3.5	A4	B1
FL Fluo Liner (1, 2, 3, 4, 5, 6)	$\frac{1}{3}$ L1 + $\frac{2}{3}$ L5	$\frac{2}{3}$ L1 + $\frac{1}{3}$ L5	1 L1	$\frac{3}{4}$ L1 + $\frac{1}{4}$ L6	$\frac{1}{2}$ L1 + $\frac{1}{2}$ L6	$\frac{1}{3}$ L2 + $\frac{2}{3}$ L5
SM Shoulder (A, B, C, D)	$\frac{1}{2}$ A + $\frac{1}{2}$ W	$\frac{2}{3}$ A + $\frac{1}{3}$ W	$\frac{3}{4}$ A + $\frac{1}{4}$ W	A 1/1	A 1/1	$\frac{1}{2}$ B + $\frac{1}{2}$ W
SM Modifier (W, T)	X	X	X	X	X	X
OD Opaque Dentine	X	X	X	X	X	X
Bleach	For Bleaching shades					
D Dentine	X	X	X	X	X	X
CD Chroma Dentine (A, B, C)	A	A	A	A	A	B
IFD Intensive Fluo Dentine (cream, yellow, orange, orange+)	Fluorescent, st					
IT Incisal transparent	IT 57	IT 57	IT 59	IT 59	IT 60	IT 57
IO Incisal opalescence (1, 2, 3, 4, 5)	Effect incisal materials for altern					
NT Neutral transparent	X	X	X	X	X	X
CM Correction						
TE Transpa Effect (1, 2, 3, 4, 5, 6)	Incisal					
G, GD Gingival	Build-up					
ST Stains universal	Modifica					

3.3. Firing test

In order to determine the temperature within your furnace, we recommend carrying out a firing test. This is the only sure way to assess your furnace's firing sequences. For the firing sample mix NT material (Neutral transparent) from Triceram® with the carving liquid LV universal.

Furnace temperatures for the firing test:

- Start temperature 500 °C / 932 °F
- Drying time 8 minutes
- Temperature rise 50 °C · min⁻¹ / 90 °F · min⁻¹
- Vacuum begin 500 °C / 932 °F
- Vacuum end once the final temperature has been met at 755 °C / 1391 °F
- Holding time 1 minute with vacuum

Place the firing sample onto platinum foil, not onto firing wool otherwise there is a risk the sample may turn cloudy. If the ceramic firing sample comes out of the furnace clear and transparent with sharp edges, then the furnace is firing correctly (Fig. 1). If the furnace is firing too high then the sample will be very shiny with no sharp edges. If it is firing too low then the sample will be milky white in colour (Fig. 2). If the furnace is firing incorrectly, increase or decrease the end temperature in 10 °C / 50 °F steps and start the test again.



Fig. 1: Optimal firing example

B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
$\frac{2}{3}$ L2 + $\frac{1}{3}$ L5	$\frac{3}{4}$ L2 + $\frac{1}{4}$ L5	1 L2	$\frac{1}{3}$ L3 + $\frac{2}{3}$ L5	$\frac{2}{3}$ L3 + $\frac{1}{3}$ L5	$\frac{3}{4}$ L3 + $\frac{1}{4}$ L5	1 L3	$\frac{2}{3}$ L4 + $\frac{1}{3}$ L5	$\frac{3}{4}$ L4 + $\frac{1}{4}$ L5	$\frac{1}{3}$ L2 + $\frac{2}{3}$ L4
$\frac{2}{3}$ B + $\frac{1}{3}$ W	B 1/1	B 1/1	$\frac{2}{3}$ C + $\frac{1}{3}$ W	$\frac{3}{4}$ C + $\frac{1}{4}$ W	$\frac{3}{4}$ C + $\frac{1}{4}$ W	C 1/1	$\frac{1}{2}$ A + $\frac{1}{3}$ D + $\frac{1}{6}$ W	$\frac{1}{2}$ A + $\frac{1}{4}$ D + $\frac{1}{4}$ W	D 1/1
X	X	X	X	X	X	X	X	X	X
X	X	X	X	X	X	X	X	X	X
not included in the Vitapan* Classic shade wheel, white and light bleaching materials									
X	X	X	X	X	X	X	X	X	X
B	B	B	C	C	C	C			
firing dentine materials to regulate the brightness within the ceramic build-up									
IT 59	IT 59	IT 59	IT 59	IT 59	IT 59	IT 60	IT 59	IT 59	IT 59
alternating layer build-up, bluish effects in reflecting light and orange coloured effects in passing light									
X	X	X	X	X	X	X	X	X	X
For small corrections (contact points, occlusion)									
alternating layer build-up, for use in incisal area, natural light refraction									
and reconstruction of gingiva and papillae, correction of gingival defects									
options for mixing with SM, OD, D, IT, IO, TE, G and stains for glazing									

* VITA is a trademark of the VITA Zahnfabrik, Bad Säckingen, Germany



Fig. 2: Firing sample after furnace temperature too low

3.4. Cleaning the furnace

Ceramic firing furnaces should be cleaned regularly to remove contamination from the inner walls of the firing chamber. For this reason we recommend:

- Regular furnace cleaning using carbon fibre platelets (REF 260-317-00)
- Base temperature: 600°C/1112°F
- Drying time: 1 minute
- Heating up speed: 100–120°C·min⁻¹/180–216°F·min⁻¹
- End temperature: 1050°C/1922°F
- Clean the firing tray including pins
- Holding time: 10 minutes

Select a firing program without vacuum.

Observe the furnace manufacturers instructions for use!

Advice:

Keep the furnace closed as much as possible in order to prevent moisture from entering the firing chamber. Always close the furnace after use, if necessary switch to the night mode.

3.5. Firing table standard program

	Base temp. *	Drying time	Heat rate	Vacuum start	Vacuum end	Final temp.	Holding time	Cooling time
Liner bake 1 + 2	500 °C 932 °F	4 min	65 °C·min ⁻¹ 117 °F·min ⁻¹	500 °C 932 °F	800 °C 1472 °F	800 °C 1472 °F	1 min under vacuum	0 min
Shoulder bake	500 °C 932 °F	6 min	55 °C·min ⁻¹ 99 °F·min ⁻¹	500 °C 932 °F	790 °C 1454 °F	790 °C 1454 °F	1 min under vacuum	0 min
Dentine bake 1 (with LV liquid universal)**	500 °C 932 °F	6 min	55 °C·min ⁻¹ 99 °F·min ⁻¹	500 °C 932 °F	760 °C 1400 °F	760 °C 1400 °F	1.5-2 min under vacuum	0 min
Dentine bake 2 (with LV liquid universal)**	500 °C 932 °F	4 min	55 °C·min ⁻¹ 99 °F·min ⁻¹	500 °C 932 °F	760 °C 1400 °F	760 °C 1400 °F	1.5-2 min under vacuum	0 min
Glaze bake***	500 °C 932 °F	2 min	55 °C·min ⁻¹ 99 °F·min ⁻¹	–	–	760 °C 1400 °F	1 min under vacuum	0 min

For the modelling liquid MV universal

	Base temp.	Drying time	Holding time at base temp.	Heat rate	Vacuum start	Vacuum end	Final temp.	Holding time	Cooling time
Dentine bake 1 (with MV liquid universal)**	500 °C 932 °F	4 min	4 min	55 °C·min ⁻¹ 99 °F·min ⁻¹	500 °C 932 °F	760 °C 1400 °F	760 °C 1400 °F	1.5-2 min under vacuum	0 min
Dentine bake 2 (with MV liquid universal)**	500 °C 932 °F	3 min	3 min	55 °C·min ⁻¹ 99 °F·min ⁻¹	500 °C 932 °F	760 °C 1400 °F	760 °C 1400 °F	1.5-2 min under vacuum	0 min

Warning:

When using the modelling liquid MV universal please ensure to follow the correct firing program. Caution, this is different to the firing instructions for modelling liquid LV universal (see above). LV and MV universal modelling liquids must not be mixed together - danger of bubbling!

The program specifications stated above are based upon fine silver calibrated furnaces. For better results select longer drying times.

Advice:

This information is to be used as an indication which must be adjusted individually to suit individual furnace manufacturers and differences in age. The program tables are based upon fine silver calibrated furnaces. All details have been carefully researched, however no responsibility is taken for the correctness of this information. Our dental technical support team is available should you have any questions regarding the firing program for your own furnace.

Hotline Tel. No. +49 72 31 / 803 - 440

3.6. Firing tables for selected furnaces

Austromat D2											
	Start	□	↑	↑	VAC	↗	END	→	(V)	↘	2
	°C/°F		min	min	%	°C·min ⁻¹ °F·min ⁻¹	°C/°F	min:sec	min:sec	1 min	2 min
Liner bake	500/ 932	1	3	0	100	65/ 117	800/ 1472	1:00	1:00	0	0
Shoulder bake	500/ 932	1	6	0	100	55/ 99	790/ 1454	1:00	1:00	0	0
Dentine 1 (with LV liquid universal)**	500/ 932	1	6	0	100	55/ 99	760/ 1400	1:30-2:00	1:30-2:00	0	0
Dentine 2 (with LV liquid universal)**	500/ 932	1	4	0	100	55/ 99	760/ 1400	1:30-2:00	1:30-2:00	0	0
Glaze bake***	500/ 932	1	2	0	0	55/ 99	760/ 1400	1:00	1:00	0	0

Austromat D4

T					01:00	min
S					03:00	min
V	500	°C			--	min
Temp. 1	800	°C	65	°C/min	01:00	min
Temp. 2	–	°C	–	°C/min	--	min
Temp. 3	–	°C	–	°C/min	--	min
VAC	800	°C	100	%	01:00	min
Liner bake						

T					01:00	min
S					06:00	min
V	550	°C			00:00	min
Temp 1	790	°C	55	°C/min	01:00	min
Temp 2	–	°C	–	°C/min	--	min
Temp 3	–	°C	–	°C/min	--	min
VAC	790	°C	100	%	01:00	min
Shoulder bake						

T					01:00	min
S					06:00	min
V	500	°C			00:00	min
Temp. 1	760	°C	55	°C/min	01:45	min
Temp. 2	–	°C	–	°C/min	--	min
Temp. 3	–	°C	–	°C/min	--	min
VAC	760	°C	100	%	01:45	min
Dentine 1 (with LV liquid universal)**						

T					01:00	min
S					04:00	min
V	500	°C			00:00	min
Temp 1	760	°C	55	°C/min	01:45	min
Temp 2	–	°C	–	°C/min	--	min
Temp 3	–	°C	–	°C/min	--	min
VAC	760	°C	100	%	01:45	min
Dentine 2 (with LV liquid universal)**						

T					01:00	min
S					02:00	min
V	500	°C			--	min
Temp. 1	760	°C	55	°C/min	01:00	min
Temp. 2	–	°C	–	°C/min	--	min
Temp. 3	–	°C	–	°C/min	--	min
VAC	760	°C	–	%	01:00	min
Glaze bake***						

* Positioning of objects only at correct standby temperatures.

** When using modelling liquid MV universal, observe different bake sequence.

*** The glaze bake can be carried out either with or without vacuum. The glaze is intensified by increasing the holding time.

Austromat M											
	START	□	↑ min	→ min	VAC LEVEL	↗	END	→ min:sec	(V) min:sec	↘ 1 min	2 ↘ min
Liner bake 1 + 2	500 °C 932 °F	1	3	0	9	65 °C·min ⁻¹ 117 °F·min ⁻¹	800 °C 1472 °F	1:00 under vacuum	d	0	0
Shoulder bake	500 °C 932 °F	1	6	0	9	55 °C·min ⁻¹ 99 °F·min ⁻¹	790 °C 1454 °F	1:00 under vacuum	d	0	0
Dentine bake 1 (with LV liquid universal)**	500 °C 932 °F	1	6	0	9	55 °C·min ⁻¹ 99 °F·min ⁻¹	760 °C 1400 °F	1:30-2:00 under vacuum	d	0	0
Dentine bake 2 (with LV liquid universal)**	500 °C 932 °F	1	4	0	9	55 °C·min ⁻¹ 99 °F·min ⁻¹	760 °C 1400 °F	1:30-2:00 under vacuum	d	0	0
Glaze bake***	500 °C 932 °F	0	2	0	***	55 °C·min ⁻¹ 99 °F·min ⁻¹	760 °C 1400 °F	1:00***	***	0	0

Austromat 3001	
Liner bake 1 + 2	C500 T60 T180 • L9 V9 T065 • C800 T60 V0 CO LO T2 C500
Shoulder bake	C500 T60 T360 • L9 V9 T065 • C790 T60 V0 CO LO T2 C500
Dentine bake 1 (with LV liquid universal)**	C500 T60 T360 • L9 V9 T055 • C760 T90-120 V0 CO LO T2 C500
Dentine bake 2 (with LV liquid universal)**	C500 T60 T240 • L9 V9 T055 • C760 T90-120 V0 CO LO T2 C500
Glaze bake ***	C500 T60 T120 • L9 T055 • C760 T60 CO LO T2 C500
The holding time for the dentine and correction bake must be between 90 and 120 seconds.	

** When using modelling liquid MV universal, observe different bake sequence.

*** The glaze bake can be carried out either with or without vacuum. The glaze is intensified by increasing the holding time.

Multimat MCI (Mach1/Mach2)

	Pre-heating temp.	Drying	Pre-heating	Vacuum	Firing time	Firing temp.	Heat rate	Vacuum
Liner bake 1 + 2	500 °C 932 °F	3 min	1 min	1 min	2 min	820 °C 1508 °F	65 °C·min ⁻¹ 117 °F·min ⁻¹	50 hPa
Shoulder bake	500 °C 932 °F	4 min	2 min	1 min	2 min	810 °C 1490 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	50 hPa
Dentine bake 1 (with LV liquid universal)**	500 °C 932 °F	4 min	2 min	1 min	2.5-3 min	780 °C 1436 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	50 hPa
Dentine bake 2 (with LV liquid universal)**	500 °C 932 °F	4 min	2 min	1 min	2.5-3 min	780 °C 1436 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	50 hPa
Glaze bake***	500 °C 932 °F	2 min	2 min	***	1.5 min	780 °C 1436 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	***

P90/P95

	Base temp.	Heat rate	Firing temp.	Closing time	Holding time	Vacuum on	Vacuum off
Liner bake 1 + 2	500 °C 932 °F	65 °C·min ⁻¹ 117 °F·min ⁻¹	800 °C 1472 °F	4 min	1 min under vacuum	500 °C 932 °F	800 °C 1472 °F
Shoulder bake	500 °C 932 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	790 °C 1454 °F	6 min	1 min under vacuum	500 °C 932 °F	790 °C 1454 °F
Dentine bake 1 (with LV liquid universal)**	500 °C 932 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	760 °C 1400 °F	6 min	1.5-2 min under vacuum	500 °C 932 °F	760 °C 1400 °F
Dentine bake 2 (with LV liquid universal)**	500 °C 932 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	760 °C 1400 °F	4 min	1.5-2 min under vacuum	500 °C 932 °F	760 °C 1400 °F
Glaze bake***	500 °C 932 °F	55 °C·min ⁻¹ 99 °F·min ⁻¹	760 °C 1400 °F	2 min	1 min***	***	–

Vacumat 100/200/50

	Final temp.	Pre-drying time	Heat rate time	Holding time	Vacuum	Cooling
Liner bake 1 + 2	800 °C 1472 °F	4 min	5 min	1 min under vacuum	6 min	–
Shoulder bake	790 °C 1454 °F	6 min	5 min	1 min under vacuum	6 min	–
Dentine bake 1 (with LV liquid universal)**	760 °C 1400 °F	6 min	5 min	1.5-2 min under vacuum	6.5-7 min	–
Dentine bake 2 (with LV liquid universal)**	760 °C 1400 °F	4 min	5 min	1.5-2 min under vacuum	6.5-7 min	–
Glaze bake***	760 °C 1400 °F	2 min	5 min	1 min***	***	–

** When using modelling liquid MV universal, observe different bake sequence.

*** The glaze bake can be carried out either with or without vacuum. The glaze is intensified by increasing the holding time.

4. Individual products

FL Fluo Liner

Fluo Liner 1	292-514-40	20 g
Fluo Liner 2	292-524-40	20 g
Fluo Liner 3	292-534-40	20 g
Fluo Liner 4	292-544-40	20 g
Fluo Liner 5	292-554-40	20 g
Fluo Liner 6	292-564-40	20 g

O Opaque

Opaque A1	292-211-30	15 g
Opaque A1	292-211-60	40 g
Opaque A2	292-212-30	15 g
Opaque A2	292-212-60	40 g
Opaque A3	292-213-30	15 g
Opaque A3	292-213-60	40 g
Opaque A3.5	292-214-30	15 g
Opaque A3.5	292-214-60	40 g
Opaque A4	292-215-30	15 g
Opaque A4	292-215-60	40 g
Opaque B1	292-221-30	15 g
Opaque B1	292-221-60	40 g
Opaque B2	292-222-30	15 g
Opaque B2	292-222-60	40 g
Opaque B3	292-223-30	15 g
Opaque B3	292-223-60	40 g
Opaque B4	292-224-30	15 g
Opaque B4	292-224-60	40 g
Opaque C1	292-231-30	15 g
Opaque C1	292-231-60	40 g
Opaque C2	292-232-30	15 g
Opaque C2	292-232-60	40 g
Opaque C3	292-233-30	15 g
Opaque C3	292-233-60	40 g
Opaque C4	292-234-30	15 g
Opaque C4	292-234-60	40 g
Opaque D2	292-242-30	15 g
Opaque D2	292-242-60	40 g

Opaque Gingiva

Opaque D3	292-243-30	15 g
Opaque D3	292-243-60	40 g
Opaque D4	292-244-30	15 g
Opaque D4	292-244-60	40 g
Opaque Gingiva	292-250-30	15 g
Chroma Opaque CO 1 Bleach	292-301-30	15 g

SM Shoulder material

Shoulder material A	293-110-30	15 g
Shoulder material B	293-120-30	15 g
Shoulder material C	293-130-30	15 g
Shoulder material D	293-140-30	15 g

SM Modifier

SM Modifier white	293-150-30	15 g
SM Modifier transparent	293-199-30	15 g

OD Opaque Dentine

Opaque Dentine A1	294-111-30	15 g
Opaque Dentine A1	294-111-60	40 g
Opaque Dentine A2	294-112-30	15 g
Opaque Dentine A2	294-112-60	40 g
Opaque Dentine A3	294-113-30	15 g
Opaque Dentine A3	294-113-60	40 g
Opaque Dentine A3.5	294-114-30	15 g
Opaque Dentine A3.5	294-114-60	40 g
Opaque Dentine A4	294-115-30	15 g
Opaque Dentine A4	294-115-60	40 g
Opaque Dentine B1	294-121-30	15 g
Opaque Dentine B1	294-121-60	40 g
Opaque Dentine B2	294-122-30	15 g
Opaque Dentine B2	294-122-60	40 g
Opaque Dentine B3	294-123-30	15 g
Opaque Dentine B3	294-123-60	40 g
Opaque Dentine B4	294-124-30	15 g
Opaque Dentine B4	294-124-60	40 g
Opaque Dentine C1	294-131-30	15 g
Opaque Dentine C1	294-131-60	40 g

OD Opaque Dentine

Opaque Dentine C2	294-132-30	15 g
Opaque Dentine C2	294-132-60	40 g
Opaque Dentine C3	294-133-30	15 g
Opaque Dentine C3	294-133-60	40 g
Opaque Dentine C4	294-134-30	15 g
Opaque Dentine C4	294-134-60	40 g
Opaque Dentine D2	294-142-30	15 g
Opaque Dentine D2	294-142-60	40 g
Opaque Dentine D3	294-143-30	15 g
Opaque Dentine D3	294-143-60	40 g
Opaque Dentine D4	294-144-30	15 g
Opaque Dentine D4	294-144-60	40 g

D Dentine

Dentine A1	294-211-30	15 g
Dentine A1	294-211-60	40 g
Dentine A2	294-212-30	15 g
Dentine A2	294-212-60	40 g
Dentine A3	294-213-30	15 g
Dentine A3	294-213-60	40 g
Dentine A3.5	294-214-30	15 g
Dentine A3.5	294-214-60	40 g
Dentine A4	294-215-30	15 g
Dentine A4	294-215-60	40 g
Dentine B1	294-221-30	15 g
Dentine B1	294-221-60	40 g
Dentine B2	294-222-30	15 g
Dentine B2	294-222-60	40 g
Dentine B3	294-223-30	15 g
Dentine B3	294-223-60	40 g
Dentine B4	294-224-30	15 g
Dentine B4	294-224-60	40 g
Dentine C1	294-231-30	15 g
Dentine C1	294-231-60	40 g
Dentine C2	294-232-30	15 g
Dentine C2	294-232-60	40 g
Dentine C3	294-233-30	15 g
Dentine C3	294-233-60	40 g
Dentine C4	294-234-30	15 g
Dentine C4	294-234-60	40 g
Dentine D2	294-242-30	15 g
Dentine D2	294-242-60	40 g
Dentine D3	294-243-30	15 g
Dentine D3	294-243-60	40 g
Dentine D4	294-244-30	15 g
Dentine D4	294-244-60	40 g

Bleach

Chroma Opaque CO 1 Bleach	292-301-30	15 g
Chroma Dentine CD 1 Bleach	294-251-30	15 g
Chroma Dentine CD 2 Bleach	294-252-30	15 g
Chroma Incisal CI 1 Bleach	295-701-30	15 g

CD Chroma Dentine

Chroma Dentine A	294-410-30	15 g
Chroma Dentine B	294-420-30	15 g
Chroma Dentine C	294-430-30	15 g

IFD Intensive Fluo Dentine

Intensive Fluo Dentine creme	294-751-30	15 g
Intensive Fluo Dentine yellow	294-760-30	15 g
Intensive Fluo Dentine orange	294-770-30	15 g
Intensive Fluo Dentine orange plus	294-773-30	15 g

IT Incisal transparent

Incisal 57	295-150-30	15 g
Incisal 57	295-150-60	40 g
Incisal 58	295-154-30	15 g
Incisal 58	295-154-60	40 g
Incisal 59	295-160-30	15 g
Incisal 59	295-160-60	40 g
Incisal 60	295-170-30	15 g
Incisal 60	295-170-60	40 g

IO Incisal opalescence

Incisal opal.1 ice white	295-210-30	15 g
Incisal opal.2 polar white	295-215-30	15 g
Incisal opal.3 snow white	295-250-30	15 g
Incisal opal.4 sky blue	295-270-30	15 g
Incisal opal.5 amber	295-280-30	15 g
Incisal opalescence 57	295-300-30	15 g
Incisal opalescence 58	295-305-30	15 g
Incisal opalescence 59	295-310-30	15 g
Incisal opalescence 60	295-315-30	15 g

NT Neutral transparent

Neutral transparent	295-299-30	15 g
Neutral transparent	295-299-60	40 g

TE Transpa Effect

Transpa Effect 1 white	295-350-30	15 g
Transpa Effect 2 yellow	295-360-30	15 g
Transpa Effect 3 pink	295-385-30	15 g
Transpa Effect 4 homey	295-370-30	15 g
Transpa Effect 5 blue	295-380-30	15 g
Transpa Effect 6 brown	295-367-30	15 g

ST Stains universal

Body Stains A	296-110-10	3 g
Body Stains B	296-120-10	3 g
Body Stains C	296-130-10	3 g
Stains 0 neutral	296-199-10	3 g
Stains 1 white	296-150-10	3 g
Stains 2 vanilla	296-162-10	3 g
Stains 3 yellow	296-160-10	3 g
Stains 45 orange plus	296-171-10	3 g
Stains 6 pink	296-185-10	3 g
Stains 7 dark pink	296-187-10	3 g
Stains 8 purple	296-189-10	3 g
Stains 9 blue	296-180-10	3 g
Stains 10 grey	296-155-10	3 g
Stains 11 olive green	296-193-10	3 g
Stains 12 olive yellow	296-192-10	3 g
Stains 13 medium brown	296-167-10	3 g
Stains 14 reddish brown	296-168-10	3 g
Stains 15 black	296-198-10	3 g

CM Correction

Correction material	297-199-10	3 g
---------------------	-------------------	-----

G Gingiva

Gingiva G (light)	298-110-30	15 g
Gingiva GD (dark)	298-120-30	15 g

Liquids


Stains universal	299-110-40	20ml
Shoulder material univ.	299-120-40	20ml
Modelling MV univ.	299-150-40	20ml
Modelling MV univ.	299-150-80	100ml
Modelling liquid LV standard univ.	299-160-41	20ml
Modelling liquid LV standard univ.	299-160-81	100ml
Modelling liquid LV standard univ.	299-160-91	500ml
Modelling liquid LV+ univ.	299-160-40	20ml
Modelling liquid LV+ univ.	299-160-80	100ml
Modelling liquid LV+ univ.	299-160-99	500ml
Universal paste liquid	299-170-40	20ml
Universal paste liquid	299-170-80	100ml
Universal B.O.L. liquid	299-180-40	20ml
Universal B.O.L. liquid	299-180-80	100ml
SM-Isofit	260-324-01	20ml
AESTHETIC Bonder liquid	299-200-10	5ml

Triceram® Starter set

290-113-02	1 piece
-------------------	---------

Triceram® Master set Zr (3)

290-112-03	1 piece
-------------------	---------



5. Description of symbols used by Dentaureum



See instructions for use



→ www.dentaureum.de

Instructions for use available in the internet

α

Symbol for the coefficient of thermal expansion CTE. If not otherwise stated, for the temperature interval from 25–500 °C / 77–932 °F after two dentine bakes.

T_g

Transformation temperature



LOT number (CH.-B)



use by ... (expiry date)

CE 0483

CE Mark MP Class 2

Rx only

Caution: Fed. Law restricts this device to sale by or on the order of a certified dental technician. (To be used exclusively by trained personnel)



Additional information regarding Dentaaurum products can be found on the internet.

www.dentaaurum.de
www.dentaaurum.de

CE 0483

Date of information: 02/08

Subject to modifications

D
DENTAURUM