




# ORIGIN<sup>®</sup>

MEETS

# VITA<sup>®</sup>

Precision Digital Solutions



Device & method to hold multiple ceramic blocks & milling toolpath are U.S. & Int'l patents pending



## With VITABLOCS<sup>®</sup>, you can't make it any quicker and easier.

Vitablocs are made from Vita's Mark II feldspathic porcelain material, which was developed in conjunction with MIT. It is the only porcelain material available that is specifically designed for highly reliable performance in CAD/CAM applications. Other materials are derivatives of conventional materials and are more susceptible to chipping and micro-cracking during the milling process. Also, Vitablocs require no crystallization firing or expensive furnace, and with the Vita Triluxe forte blocs there is no need for staining and glazing in order to get a highly esthetic result. Users simply mill to full contour, polish, and place the restoration in the mouth. Simplicity, time and cost savings - the foundation of the ORIGIN systems.

More than 25 years of clinical experience

More than 16 million restorations

Unsurpassed clinical survival rates\*


\*85% at 18 years; as "Good as gold", 94% at 7 years for VITA Mark II - Data on File.

# VITABLOCS<sup>®</sup>



### Cutting Titanium!

- ORIGIN NobleTitan™ Titanium
- Titanium is one of the most bio-compatible metals available.
- ADA classifies titanium between high noble and noble alloys.
- Ideal CAD/CAM substitute for base metal and noble metal.
- Titanium porcelain available.



## Most Cost Efficient Production

Premium OPEN Source Dental CAD/CAM Systems at Competitive "OPEN MARKET" Prices

### Zirconia Sintering in 2 Hours

Non-microwave Breakthrough Technology!  
U.S. & International Patent Pending

### Comprehensive System Evaluation & Training

Call us to schedule a system demonstration and receive hands-on training at our facility in Salt Lake City, Utah. Bring your own cases to scan, design, mill and sinter. Spend as much time as you need to thoroughly evaluate and test the ORIGIN CAD/CAM system, we will cover your meals and hotel stay.

### Best of Class Components

Each stand-alone ORIGIN component (scanner, milling machine, sintering furnace, YZ zirconia & software) is among the best in its field. Each component has a proven track record of performance and dependability, allowing you to trust that your system will provide the exacting results that you demand.

### Korean Technicians

### 한국어 상담 및 서비스

PFM accounts들을 Zirconia crown 으로 교체하기에 지금보다 더 좋은 시기는 없습니다. 장비에 관한 상담 및 In-house milling 에 관한 know-how 등 일체의 내용을 한국어로편하게 서비스 해 드립니다. 업계최초의 자신있는 coping 원가공개 및 다른 system 과는 월등하게 차별화된 서비스로 제공되는 ORIGIN CAD CAM 대열에 참여하십시오.

(DuoTron Oven은 한국의 신기술로 제작 되었으며 나머지 제품들은 모두 검증받은 유럽의 Best Seller제품입니다.)

www.OriginCadCam.com  
support@OriginCadCam.com

800-255-2839

2371 So. Presidents Dr., Ste. E & F  
West Valley City, UT 84120

